ACTH

¹²⁵I RIA Kit

For the quantitative determination
of human ACTH in EDTA plasma

Instruction Manual

Manuel d'Instructions

Testanleitung

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REF: 24130

DiaSorin

Stillwater, Minnesota 55082-0285, U.S.A.

TABLE OF CONTENTS

English	Page 1
Français	12
Deutsch	23
Italiano.....	35

ADRENOCORTICOTROPIC HORMONE RADIOIMMUNOASSAY

1. INTENDED USE

FOR *IN VITRO* DIAGNOSTIC USE.

This kit contains instructions and materials for the quantitative determination of human ACTH in EDTA plasma by radioimmunoassay (RIA).

2. SUMMARY AND EXPLANATION

ACTH is a 39 amino acid pituitary hormone with the following structure:

NH₂-Ser-Tyr-Ser-Met-Glu-His-Phe-Arg-Trp-Gly-Lys-Pro-Val-
Gly-Lys-Lys-Arg-Arg-Pro-Val-Lys-Val-Tyr-Pro-Asn-Gly-
Ala-Glu-Asp-Glu-Ser-Ala-Glu-Ala-Phe-Pro-Leu-Glu-Phe-OH

The sequence 1-13 is the same as alpha melanocyte stimulating hormone. Sequence 1-24 is considered the active portion of the molecule. Sequence 18-39 is known as clip or corticotropin-like intermediate lobe peptide.¹ The molecule is generally a basic protein and quite labile to the action of peptidases. Since ACTH readily binds to glass and associates with serum proteins, a special buffer preparation is necessary to ensure preservation of ACTH and to reduce nonspecific binding.^{1,2}

ACTH stimulates the secretion of cortisol from the adrenal gland and is under control of a feedback mechanism, forming a tightly controlled loop between the hypothalamus, pituitary, and adrenal glands.^{3,4}

Two factors in the hypothalamus will stimulate ACTH release. One is vasopressin, a well known peptide, and the other is a lesser known substance, CRF or corticotropin releasing factor. Vasopressin infusion is a potentially useful provocative test; however, it has certain dangers and must be used cautiously.

ACTH measurement is a clinically useful tool in the assessment of adrenal and pituitary insufficiency.⁵ It is also useful in assessment of Cushing's syndrome, whether related to pituitary or adrenal tumors.^{6,7} Further information on the use of the DiaSorin ACTH assay is available upon request.

3. PRINCIPLE OF THE ASSAY

The DiaSorin ACTH assay is an equilibrium radioimmunoassay which can be performed by one of two methods.

In Option A, sample, first antibody and tracer are combined and incubated for 16-24 hours at 2-8°C. A pre-precipitated second antibody complex is added to separate the bound from free tracer. The assay can then be centrifuged and decanted after a 15 minute incubation at 20-25°C.

In Option B, sample, first antibody, and tracer are combined and incubated for 3 hours at 2-8°C. A pre-precipitated second antibody complex is added to separate the bound from free tracer. The assay can then be centrifuged and decanted after a 15 minute incubation at 20-25°C. The 3 hour option is less sensitive than the overnight procedure, and the results obtained will be slightly higher.

Quality control will be performed at DiaSorin using only Option A.

4. REAGENT PROVIDED IN THE KIT

ACTH Calibrator 0	1 vial/3 mL	1 vial/3 mL
ACTH Calibrators (1-5)	5 vials/1.0 mL	5 vials/1.0 mL
ACTH Antiserum	1 vial/14 mL	2 vials/14 mL
¹²⁵ I ACTH	1 vial/14 mL	2 vials/14 mL
Precipitating Complex (GAR-PPT)	1 vial/ 35 mL	2 vials/ 35 mL
ACTH Controls level 1 and 2	2 vials/1.0 mL	2 vials/1.0 mL
Number of tests	65	130

STORAGE: Upon receipt, and prior to reconstitution, store all reagents at 2-8°C. After reconstitution store all reagents at -15° or lower until the expiration date on the label. Reagents should not be used past the expiration date.

When reconstituting the contents of the vials, mix gently to avoid foaming. Reagents from different batches must not be mixed.

4.1 Calibrator 0: lyophilized reagent

ACTH-free human plasma with 0.02% sodium azide and other stabilizers added. Reconstitute the vial with 3.0 mL of purified water and let stand on crushed ice for 15-20 minutes until completely dissolved; mix thoroughly before using.

4.2 ACTH Calibrators (A-E): lyophilized reagent

ACTH calibrators (porcine), at nominal concentrations ranging from 20-500 pg/mL (4.5-112 picomoles/L), are pre-diluted in Calibrator 0 containing 0.02% sodium azide and other stabilizers. Exact concentration values are assigned with each lot. Reconstitute each vial with 1.0 mL of purified water and let stand on crushed ice for 15-20 minutes until completely dissolved; mix thoroughly before using. For consistent results, these calibrators can be thawed only once for use in an assay. The DiaSorin ACTH calibrator has been calibrated against the World Health Organization (W.H.O.) standard 74/555. Any comparison with other products or procedures should be done with this reference calibrator. The calibrator demonstrates commutability with patient samples when used with reagents and operating procedure of this *in vitro* diagnostic test as recommended.

4.3 ACTH Antiserum: lyophilized reagent

Rabbit anti-ACTH serum diluted in BSA-borate buffer with 0.002% Thimerosal and other stabilizers added. Reconstitute each vial with 14 mL of purified water and let stand on crushed ice for 15-20 minutes until completely dissolved; mix thoroughly before using.

4.4 ¹²⁵I ACTH: lyophilized reagent

Synthetic ACTH labeled with iodine-125, diluted in BSA-borate-EDTA buffer with 0.008% Thimerosal and other stabilizers added. Reconstitute each vial with 14 mL of purified water and let stand on crushed ice for 15-20 minutes until completely dissolved; mix thoroughly before using. For consistent results, thaw only once for use in an assay.

4.5 Precipitating Complex (GAR-PPT): lyophilized reagent

Normal rabbit serum, pre-precipitated with goat anti-rabbit serum and polyethylene glycol (PEG), is diluted in BSA-borate buffer with 0.03% Thimerosal added. Reconstitute the vial with 35 mL of purified water. Mix THOROUGHLY until the suspension appears homogeneous and then allow the vial to stand for a minimum of 30 minutes at room temperature with occasional mixing.

4.6 ACTH Control (Level 1 and 2): lyophilized reagent

Human plasma spiked, if necessary, with the appropriate amount of synthetic ACTH to obtain a concentration within a specified range. 0.1% Sodium azide and other stabilizers are added. Reconstitute each vial with 1.0 mL of purified water and place on crushed ice for 15-20 minutes until completely dissolved; mix thoroughly and treat as an unknown sample. For consistent results this control can be thawed only once for use in an assay. Values for each lot are printed on the control vial.

5. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE.

Not for internal or external use in humans or animals.

REAGENTS CONTAINING HUMAN SOURCE MATERIAL

Treat as potentially infectious.

Each serum/plasma donor unit used in the preparation of this product has been tested by a U.S. FDA approved method and found non-reactive for the presence of HBsAg, antibody to HCV, and antibody to HIV 1/2. While these methods are highly accurate, they do not guarantee that all infected units will be detected. This product may also contain other human source material for which there is no approved test. Because no known test method can offer complete assurance that hepatitis B virus (HBV), hepatitis C virus (HCV), Human Immunodeficiency Virus (HIV), or other infectious agents are absent, all products containing human source material should be handled in accordance with good laboratory practices using appropriate precautions as described in the U.S. Centers for Disease Control and Prevention/National Institutes of Health Manual, "Biosafety in Microbiological and Biomedical Laboratories," 4th ed., May 1999 or current edition.

REAGENTS CONTAINING SODIUM AZIDE

CAUTION: Some reagents in this kit contain sodium azide. Sodium azide may react with lead or copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up. For further information, refer to "Decontamination of Laboratory Sink Drains to Remove Azide Salts," in the Manual Guide-Safety Management No. CDC-22 issued by the Centers for Disease Control and Prevention, Atlanta, GA, 1976.

European Communities Hazardous Substance Risk Phrases (Council Directive 1999/45/EC)

R20/21/22 - Harmful by inhalation, in contact with skin and if swallowed.

R32 - Contact with acids liberates very toxic gas.

S28 - After contact with skin, wash immediately with plenty of water.

REAGENTS CONTAINING THIMEROSAL

Some reagents in this kit contain thimerosal which contains a mercury compound. Disposal of elemental mercury, inorganic mercury, mercury oxides, and mercury compounds should be done in strict compliance with all local, state, and federal regulations.

WARNING: This product contains a chemical known to the State of California to cause birth defects or other reproductive harm.

REAGENTS CONTAINING IODINE-125

This kit contains radioactive material which does not exceed 1.5 μCi (55.5 kBq) for the 65 tube kit and 3 μCi (111 kBq) for the 130 tube kit of iodine-125. Appropriate precautions and good laboratory practices should be used in the storage, handling, and disposal of this material.

For practitioners or institutions receiving radioisotopes under a general license:

This radioactive material may be received, acquired, possessed, and used only by physicians, veterinarians in the practice of veterinary medicine, clinical laboratories or hospitals, and only for *in vitro* clinical or laboratory tests not involving internal or external administration of the material, or the radiation therefrom, to human beings or animals. Its receipt, acquisition, possession, use, and transfer are subject to the regulations and the general license of the U.S. Nuclear Regulatory Commission or of the state with which the Commission has entered into an agreement for the exercise of regulatory authority.

1. Storage of radioactive material should be limited to a specifically designated area.
2. Access to radioactive materials must be limited to authorized personnel only.
3. Do not pipette radioactive material by mouth.

4. Do not eat or drink within designated radioactive work areas.
5. Areas where spills may occur should be wiped up, then washed with an alkali detergent or radiological decontamination solution. Any glassware used must be rinsed completely with water before washing with other laboratory glassware.

For practitioners or institutions receiving radioisotopes under a specific license:

The receipt, use, transfer, and disposal of radioactive materials are subject to the regulations and conditions of your specific license.

WARNING: This product contains a chemical known to the State of California to cause cancer.

ATTENTION: Radioactivity printed in the package insert may be slightly different from the radioactivity printed on the box label and on the tracer vial label. The box label and the tracer vial label indicate the actual amount of radioactivity at the calibration date where the package insert indicates the theoretical radioactivity of the kit.

6. INDICATIONS OF POSSIBLE DETERIORATION OF KIT REAGENTS

- 6.1 The presence of abnormal particulate matter in any of the reagents.
- 6.2 A shift in the slope or position of the calibrator curve from what is normally obtained.
- 6.3 A decrease in maximum binding.
- 6.4 A high nonspecific binding.

7. COLLECTION AND PREPARATION OF SPECIMEN

One hundred microliters in duplicate of EDTA plasma are required for the assay.

Collect blood by venipuncture in a 5 or 10 mL evacuated glass tube. EDTA (7.2 mg/ 5 mL whole blood) is used as anticoagulant. Mix gently and place on ice. Centrifuge for 15 minutes at $760 \times g$ to obtain hemolysis-free plasma. A short centrifugation at 20°C will not be detrimental; however, refrigerated centrifugation is recommended. All plastics, glassware, or other material coming into contact with the specimen should be entirely free of any contamination.

PRECAUTION: Heparinized plasma yields falsely low values.

Enzyme inhibitors are unnecessary at the sample collection phase; however, it is important to keep the sample cold ($0-4^\circ\text{C}$) after collection (for up to 8 hours). For long-term storage add 500 KIU/mL of Aprotinin to the specimen and store it at -15°C . It is very important to keep the specimen frozen, preferably in plastic tubes. A decrease in ACTH may be experienced with repeated freeze-thaw of the specimen.

8. EQUIPMENT AND MATERIALS REQUIRED, BUT NOT SUPPLIED

- 8.1 Disposable borosilicate glass tubes, 12 x 75 mm.
- 8.2 Temperature controlled centrifuge to accommodate 12 x 75 mm tubes ($20-25^\circ\text{C}$).
- 8.3 Gamma scintillation counter capable of counting ^{125}I .
- 8.4 Vortex.
- 8.5 Pipetting devices:
 - a. Micropipettors calibrated to deliver 50 μL , 100 μL and 200 μL .
 - b. Repeating dispensers, calibrated to deliver 200 μL and 500 μL .
- 8.6 Purified water.

9. ASSAY PROCEDURE - OPTION A (OVERNIGHT)

- 9.1 Thaw unknown samples and place on crushed ice.
- 9.2 Reconstitute the lyophilized reagents and allow any frozen reagents to thaw completely. Do not allow reagents to reach temperatures above 25°C . Mix reagents gently and place on crushed ice before using.

$$* g = (1118 \times 10^{-8}) (\text{radius in cm}) (\text{rpm})^2$$

- 9.3 Set up labeled 12 x 75 mm disposable glass tubes in duplicate according to the Scheme of the Assay, on the back page.
 - 9.4 Place the rack of tubes on crushed ice.
 - 9.5 Add reagents as follows:
 - a. **Total count tubes**
Set aside until Step 6
 - b. **Nonspecific binding tubes (NSB)**
100 μ L of Calibrator 0
 - c. **Calibrator 0**
100 μ L of Calibrator 0
200 μ L of ACTH antiserum
 - d. **ACTH Calibrators (A-E)**
100 μ L of ACTH calibrator
200 μ L of ACTH antiserum
 - e. **Controls plasma and unknown samples**
100 μ L of plasma
200 μ L of ACTH antiserum
 - 9.6 Add 200 μ L of 125 I ACTH to all tubes.
 - 9.7 Mix and incubate for 16-24 hours at 2-8°C.
 - 9.8 Vigorously mix the GAR-PPT; add 500 μ L to all the tubes except the total count tubes.
 - 9.9 Mix and incubate for 15-25 minutes at 20-25°C.
 - 9.10 Centrifuge the tubes for 20 minutes using 760 x g* at 20-25°C.
 - 9.11 Immediately decant the supernatant from all the tubes except the total count tubes by inverting them for a minimum of 2 minutes. Blot the tubes on absorbent paper to remove any drops of supernatant that may be remaining on the rims before turning the tubes upright.
 - 9.12 Using a gamma scintillation counter, count the precipitate of each tube and the total count tubes for a sufficient time to achieve statistical accuracy. (See Limitations of the Procedure.)
- 10. ASSAY PROCEDURE - OPTION B (3 HOUR ASSAY)**
- 10.1 Thaw unknown samples and place on crushed ice.
 - 10.2 Reconstitute the lyophilized reagents and allow any frozen reagents to thaw completely. Do not allow reagents to reach temperatures above 25°C. Mix reagents gently and place on crushed ice before using.
 - 10.3 Set up labeled 12 x 75 mm disposable glass tubes in duplicate according to the Scheme of the Assay, on the back page.
 - 10.4 Place the rack of tubes on crushed ice.
 - 10.5 Add reagents as follows:
 - a. **Total count tubes**
Set aside until Step 6
 - b. **Nonspecific binding tubes (NSB)**
100 μ L of Calibrator 0
 - c. **Calibrator 0**
100 μ L of Calibrator 0
200 μ L of ACTH antiserum
 - d. **ACTH Calibrators (A-E)**
100 μ L of ACTH calibrator
200 μ L of ACTH antiserum

$$* g = (1118 \times 10^{-8}) (\text{radius in cm}) (\text{rpm})^2$$

- e. **Controls plasma and unknown samples**
 - 100 µL of plasma
 - 200 µL of ACTH antiserum
- 10.6 Add 200 µL of ¹²⁵I ACTH to all tubes.
- 10.7 Mix and incubate for 3 hours (±15 minutes) at 2-8°C.
- 10.8 Vigorously mix the GAR-PPT; add 500 µL to all the tubes except the total count tubes.
- 10.9 Mix and incubate for 15-25 minutes at 20-25°C.
- 10.10 Centrifuge the tubes for 20 minutes using 760 x g* at 20-25°C.
- 10.11 Immediately decant the supernatant from all the tubes except the total count tubes by inverting them for a minimum of 2 minutes. Blot the tubes on absorbent paper to remove any drops of supernatant that may be remaining on the rims before turning the tubes upright.
- 10.12 Using a gamma scintillation counter, count the precipitate of each tube and the total count tubes for a sufficient time to achieve statistical accuracy. (See Limitations of the Procedure.)

11. PROCEDURAL COMMENTS

- 11.1 Add each aliquot of reagent to the lower third of the assay tube to ensure complete mixture of reagents.
- 11.2 Some manufacturers' disposable borosilicate glass tubes yield elevated nonspecific bindings.
- 11.3 If you choose to aspirate the supernatant from the precipitate, be careful not to disturb the precipitate.
- 11.4 To completely monitor the consistent performance of an RIA there are additional factors which must be checked. DiaSorin suggests a check of the following parameters to ensure consistent kit performance.
 - a. **Total Counts**
 - b. **Maximum Binding**
Average counts per minute (CPM) of Calibrator 0 Tube / Average CPM of Total Count Tubes.
 - c. **Nonspecific Binding**
Average CPM of NSB Tube / Average CPM of Total Count Tubes.
 - d. **Slope of Calibrator Curve**
For example, monitor the 80, 50, and 20% points of the calibrator line.

12. QUALITY CONTROL

Each laboratory should include control plasma in every assay to ensure the validity of each assay's results. A mean and standard deviation should then be determined for each control using a minimum of ten assays. An acceptable range of values may then be obtained for these controls using ±2 standard deviations of the values previously determined. The DiaSorin Quality Control Laboratory has determined a range for the controls included in this kit using Option A.

13. CALCULATION OF RESULTS

There are many methods in existence for calculating results of RIAs. Each is based on obtaining a calibration curve by plotting the extent of binding against stated concentrations of the calibration calibrators. This graph may be either linear or logarithmic scale. Each of these methods gives essentially the same values for controls and samples, although certain assays may "fit" better into one particular method versus another. The calculation method for DiaSorin Quality Control Laboratory is % B/B₀ versus log concentration.

$$* g = (1118 \times 10^{-8}) (\text{radius in cm}) (\text{rpm})^2$$

- 13.1 Calculate the average CPM for each calibrator, control and unknown sample.
- 13.2 Subtract the average CPM of the NSB tubes from all counts.
- 13.3 Divide the corrected CPM of each calibrator, control, or unknown sample by the corrected CPM of the Calibrator 0.

$$B/B_0 (\%) = \frac{\text{CPM of Calibrator or Unknown Sample} - \text{CPM of NSB}}{\text{CPM of Calibrator 0} - \text{CPM of NSB}} \times 100$$

- 13.4 Using 2 cycle semi-log or log-logit graph paper, plot percent B/B₀ for the ACTH calibrators (vertical axis) versus the concentration (horizontal axis).
- 13.5 Draw a best-fit line through the points.
- 13.6 Interpolate the levels of ACTH in the unknown samples from the plot.
- 13.7 If any unknown sample reads greater than the highest calibrator, it should be diluted appropriately with Calibrator 0 and assayed again.
- 13.8 If an unknown sample has been diluted, correct for the appropriate dilution factor.
- 13.9 Calculate maximum binding by dividing CPM of Calibrator 0 by the average total counts obtained in the total count tubes.

TABLE III
DiaSorin ACTH RIA Sample Data

Tube	Duplicate CPM	Average CPM	Corrected CPM	Percent Bound (B/T)	Percent (B/B ₀)	Conc. (pg/mL)
Total Count	18,327 18,750	18,539				
NSB	1,790 1,852	1,821		9.8		
Calibrator 0	9,289 9,402	9,346	7,525	50.4	100	
Calibrators (pg/mL)						
A (20)	8,765 8,721	8,743	6,922		92.0	
B (50)	7,785 8,228	8,007	6,186		82.2	
C (100)	6,766 6,696	6,831	4,910		65.2	
D (200)	5,171 5,190	5,181	3,360		44.7	
E (500)	3,296 3,416	3,356	1,535		20.4	
Unknown Samples						
1	7,253 7,362	7,308	5,487		72.9	74
2	5,860 5,885	5,873	4,052		53.8	147

Typical sample data and a calibrator curve for Option A are shown in TABLE III and FIGURE 1; this information is for reference only and should not be used for the calculation of any value.

ACTH SAMPLE CALIBRATOR CURVE


FIGURE 1

REDUCTION DATA

The DiaSorin QC lab uses a smoothed spline curve fit.

14. LIMITATIONS OF THE PROCEDURE

- 14.1 Control values have been established using Option A, the overnight assay procedure. Values obtained using the three hour procedure may give slightly different results.
- 14.2 Counting times should be sufficient to prevent statistical error (for example, accumulation of 2,000 CPM will yield 5% error; 10,000 CPM will yield 1% error).

15. EXPECTED VALUES

15.1 Normal Range

Each laboratory should establish its normal range. The normal values observed in the DiaSorin Laboratory are based on comparison of values to WHO ACTH standard 74/555 and not to any other commercial calibrator. The normal range in EDTA plasma found in the DiaSorin Laboratory was: 0-71 pg/mL for the Option A (n=50), 0-82 pg/mL for the Option B (n=49). Since our standardization was done with EDTA plasma, the use of serum or any other type of plasma may lead to a higher or lower normal range.

15.2 Nonspecificity

Since nonspecificities can occur in any radioimmunoassay, the laboratory should always be prepared to serially dilute a sample to test for parallelism or to do an affinity gel extraction of the sample to verify whether the measured substance is immunoadsorbable. The physician, in turn, should be aware that cortisol measurements, provocative and/or suppression tests are often performed with ACTH to evaluate the patient's pituitary-adrenal state.

16. SPECIFIC PERFORMANCE CHARACTERISTICS

16.1 Precision

Intra Assay Variation (values = pg/mL)

Sample Number	Mean Value	S.D	%C.V.
Low	69	8.6	12.5
Medium	281	17.6	6.3
High	600	48.8	8.1

Inter Assay Variation (values = pg/mL)

Sample Number	Mean Value	S.D	%C.V.
Low	23	1.4	6.0
High	218	14.7	6.7

16.2 TRUENESS: THE ASSAY TRUENESS HAS BEEN VERIFIED BY THE DILUTION TEST AND THE RECOVERY TEST.

Linearity (Parallelism)

Serial Dilution Study of four Patient Samples (values = pg/mL)

Sample Number	Undiluted	1/2	1/4	1/8
1	156	161	158	--
2	483	476	555	--
3	>500	---	988	982
4	>500	---	930	1008

Recovery (Accuracy) (values = pg/mL)

Background	Calibrator Added	Expected Value	Measured Value	Percent Recovery
Set No. 1 53	37.5	90.5	81	90
	150	203	182	90
	300	353	276	78
Set No. 2 60	37.5	97.5	90	92
	150	210	190	90
	300	360	315	88
Set No. 3 40	37.5	77.5	65	84
	150	190	155	82
	300	340	345	101

16.3 Analytical Sensitivity

When defined as the apparent concentration at 3 standard deviations from the counts at maximum binding, the minimum detectable amount is 15 pg/mL for the overnight option and 20 pg/mL for the 3 hour option.

16.4 Analytical Specificity

Comparison of the cross-reactivity of ACTH antibody was made with the following peptides at 1200 pg/mL:

Peptide	%Cross-reactivity
a-MSH	<0.01
b-Endorphin	<0.01
b-Lipotropin	<0.01
Leucine Enkephalin	<0.01
Methionine Enkephalin	<0.01
Bombesin	<0.01
Calcitonin	<0.01
Parathyroid Hormone	<0.01
FSH	<0.01
Human Growth Hormone	<0.01
Vasopressin	<0.01
Oxytocin	<0.01
Substance P	<0.01
Porcine ACTH 1-39	100
Human ACTH 1-24	100

Our evidence indicates that rat ACTH is also measurable with this system; however, rat plasma may contain other interfering substances.

REFER TO LAST PAGE FOR REFERENCES

SCHEME OF THE ASSAY

Option A

1. Reconstitute the lyophilized reagents and allow any frozen specimens to thaw completely. Do not allow reagents to reach above 25°C. Identify tubes in duplicate.
2. Dispense reagents according to the following scheme, on crushed ice.

Tubes/Reagents	Total Counts	NSB	Cal 0-5	Controls and unknown samples
Calibrators (0-5)	-	100 µL	100 µL	-
Controls	-	-	100 µL	-
Unknown Samples	-	-	-	100 µL
Antiserum	-	-	200 µL	200 µL
Tracer	200 µL	200 µL	200 µL	200 µL

3. Cover the tubes with parafilm and vortex gently.
4. Incubate for 16-24 hours at 2-8°C.
5. Dispense 500 µL of the GAR-PPT into all wells except the total count tubes.
6. Cover the tubes and vortex gently.
7. Incubate for 15-25 minutes at 20-25 °C.
8. Centrifuge using 760 x g* for 20 minutes.
9. Decant the supernatants.
10. Count each tube in a gamma counter for 60 seconds or longer.

Option B

1. Reconstitute the lyophilized reagents and allow any frozen specimens to thaw completely. Do not allow reagents to reach above 25°C. Identify tubes in duplicate.
2. Dispense reagents according to the following scheme, on crushed ice.

Tubes/Reagents	Total Counts	NSB	Cal 0-5	Controls and unknown samples
Calibrators (0-5)	-	100 µL	100 µL	-
Controls	-	-	100 µL	-
Unknown Samples	-	-	-	100 µL
Antiserum	-	-	200 µL	200 µL
Tracer	200 µL	200 µL	200 µL	200 µL

3. Cover the tubes with parafilm and vortex gently.
4. Incubate for 3 hours +/- 15 minutes at 2-8°C.
5. Dispense 500 µL of the GAR-PPT into all wells except the total count tubes.
6. Cover the tubes and vortex gently.
7. Incubate for 15-25 minutes at 20-25 °C.
8. Centrifuge using 760 x g* for 20 minutes.
9. Decant the supernatants.
10. Count each tube in a gamma counter for 60 seconds or longer.

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$$* g = (1118 \times 10^{-8}) (\text{raggio in cm}) (\text{rpm})^2$$

REFERENCES/BIBLOGRAPHIE/LITERATUR/BIBLIOGRAFIA


1. Yalow, R.S., "Big ACTH and Bronchogenic Carcinoma," **Annual Review of Medicine**, 30:241, (1979).
2. Orth, D.N., "Adrenocorticotrophic Hormone (ACTH)," in **Methods of Hormone Radioimmunoassay**, 2nd Edition, Academic Press, (New York), p. 245, (1979).
3. Krieger, D.T., "Factors Influencing Circadian Periodicity of ACTH and Corticosteroids," **The Medical Clinics of North America**, 62:251, (1978).
4. Krieger, D.T., "Rhythms of CRF, ACTH and Corticosteroids," **Endocrine Rhythms**, ed. D.T. Krieger, Raven Press (New York), p. 123, (1979).
5. Scott, R.S., R.A. Donald and E.A. Espiner, "Plasma ACTH and Cortisol Profiles in Addisonian Patients Receiving Conventional Substitution Therapy," **Clinical Endocrinology**, 9:571, (1978).
6. Singer, W., K. Kovacs, N. Ryan and E. Horvath, "Ectopic ACTH Syndrome; Clinico-pathological Correlations," **Journal of Clinical Pathology**, 31:591, (1978).
7. Yalow, R.S., "Ectopic ACTH in Carcinoma of the Lung," **Lung Cancer: Progress in Therapeutic Research**, ed. F.M. Muggia and M. Rozenzweig, Raven Press (New York), p. 209, (1979).

SYMBOLS USED WITH DEVICES

	English	Français	Deutsch	Italiano
	European Conformity	Conformité aux normes européennes	CE – Konformitätskennzeichnung	Conformità europea
	Expiration Date	Date limite d'utilisation	Mindesthaltbarkeitsdatum	Data di scadenza
	Manufacturer	Fabricant	Hersteller	Fabbricante
	Consult Instructions for Use	Consulter les instructions d'utilisation	Gebrauchsanweisung beachten	Consultare le istruzioni per l'uso
	Maximum temperature	Température maximale	Höchsttemperatur	Temperatura massima
IVD	In vitro diagnostic.	Diagnostic in vitro.	In-vitro-Diagnostikum.	Diagnostica in vitro.
LOT	Lot No.	No. de lot	Chargen-Nr.	Lotto n°
Ab	Antiserum	Antisérum	Antiserum	Antisiero
Ab PEG	Precipitating reagent	Réactif précipitant	Fällungsreagenz	Reagente precipitante
Ag ¹²⁵ I	Tracer: antigen labelled with ¹²⁵ I	Traceur : antigène marqué à l'iode ¹²⁵	Tracer: ¹²⁵ I-markiertes Antigen	Tracciatore: antigene etichettato con ¹²⁵ I
CAL	Calibrator	Étalon	Kalibrator	Calibratore
CONTROL	Control serum	Sérum de contrôle	Kontrollserum	Siero di controllo
	Radioactive	Radioactif	Radioaktiv	Radioattivo
	Harmful	Nocif	Gesundheitsschädlich	Nocivo

DiaSorin Inc.
1951 Northwestern Avenue
P.O. Box 285
Stillwater, MN 55082-0285

DiaSorin S.p.A.,
13040 Saluggia (VC) Italy

Phone: 651-439-9710
FAX: 651-351-5669

For Customer Service in the
U.S. and Canada Call Toll Free:
800-328-1482

In the United Kingdom Call:
+44(0) 1344 401 430
FAX: +44(0) 7884 050812

10398

35269
12/10

PRINTED IN U.S.A.